Piezo[®]

User Guide

Slide on your waist & go...

◀ Press & hold to reset your steps to zero

Setting your step goal:

What is your current average daily step count?

Track your steps for 5-7 days to find out!

Set a goal to gradually increase your daily steps by 500/day

Once you reach this goal set a new one!

Your eventual goal?

8000 steps a day

Interested in weight loss?

Research indicates you should aim for **12-15,000 steps a day**

Tips:

- To prevent loss, clip your StepStrap to a belt loop or pocket
- Change the battery when the display fades
- Your piezo uses LR-44 batteries (available in stores, or at StepsCount)

About the Piezo®

Why a Piezo?

The Piezo series activity meters are designed with a medical grade piezoelectric sensor. The Piezo has been developed for use in research and clinics. Piezos have an unmatched level of accuracy across a wide range of ages, body sizes and walking speeds.

Some interesting Piezo facts:

The Piezo is equipped with a 4 second time filter. You have to be moving for 4 continuous seconds before the steps are added.

If a Piezo does not sense any movement for 5 minutes, the LCD display will automatically shut off to conserve battery life. When you start moving again, the LCD will automatically come back on.

Tips for getting started with your Piezo:

Place your Piezo anywhere on a waistband or belt. Slide one corner of the clip onto your waist band or belt and position it. The Piezo can also be worn on a bra, on the bottom band. Regardless of where you wear it, make sure the Piezo remains horizontal (do not position on its side).

A StepStrap® is attached to every piezo. It prevents your Piezo from dropping in a puddle (or worse), should it come off your waist.

The StepStrap has saved the day more than once!

When you receive your Piezo, you can remove the protective film on the LCD. Do not immerse your pedometer in water, or apply excessive force as this may damage the LCD.

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